

WT-03 BASE

WT-04 STANDARD

WT-06 FLEX

DIGITAL THERMOMETER



EN Instructions for use

INTRODUCTION

Dear consumer! We congratulate you with a B.Well digital thermometer buying! Thank you for choosing our product!

INTENDED PURPOSE

Thermometers are used for measuring body temperature by oral, axillar and rectal method.

PRECAUTIONS


- The thermometer may be used only for measuring body temperature!
- The minimum measurement time until the beep is heard must be maintained without exception!
- Ensure that children do not use the thermometer unsupervised; some parts are small enough to be swallowed.
- Protect the thermometer from impact and dropping!
- Avoid bending the tip more than 45°!
- Avoid ambient temperatures above 60 °C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section "Cleaning and Disinfecting" to clean the thermometer.
- Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.
- Under any circumstances, the temperature taking result is ONLY for reference. Before taking any medical action, please consult your physician.

WHAT IS A "NORMAL" TEMPERATURE?

Although the generally accepted "normal" temperature reading is 37.0°C, temperature readings can vary from 36.1°C to 37.2°C and still may be considered "normal". Variations in temperature can be attributed to activities such as exercise, smoking, eating and drinking. Your temperature is lower in the morning than in the afternoon. Other variations may be due to the location of the temperature reading: Axillary (armpit measuring), Oral (mouth measuring), Rectal (rectum measuring). While the oral temperature follows the guidelines previously described, a rectal temperature is generally 0.5°C higher. Axillary temperature will be 0.5°C lower.

FEATURES OF YOUR DIGITAL THERMOMETER

WT-03 base WT-04 standard WT-06 flex



1. ON/OFF button
2. Probe tip
3. Battery cover
4. Display window

CLEANING INSTRUCTIONS

- 1st method.** Clean the thermometer with a soft cloth and isopropyl alcohol diluted with water.
- 2d method.** Wash it with warm soapy water, then wipe with a clean wet cloth and let it dry.
- 3d method.** Clean the thermometer with a soft alcohol wet wipe (wipes for injections are preferred).

PREPARING FOR USE

Press ON/OFF button to activate. The unit would beep and a display would show all the segments for around 2 seconds.

WT-03: Then the thermometer would show a test temperature which is 36,5°C.

WT-04 and WT-06: The temperature taken last time would be displayed. Then a test temperature which is 37,0°C will appear for 1 second. If the room temperature is <32,00°C, Lo°C would be displayed.

WT-03: If the room temperature is ≥42,0°C, Hi °C would be displayed.

WT-04, WT-06: If the room temperature is ≥43,0°C, Hi °C would be displayed. When °C appears, the thermometer is ready for temperature measuring.

WT-03: For viewing the temperature taken last time press and hold ON/OFF button for 2 sec.

WT-04: Press "ON/OFF" button over 2 seconds to switch between °C/°F. The display will show set unit °C or °F until the button is released. Afterwards the display will show Lo °C or Lo °F

HOW TO TAKE YOUR TEMPERATURE

• Taking an axillary temperature.

Wipe the underarm with a dry towel and place the probe tip under the arm.

⚠ **Attention!** Severe sweating may cause incorrect measuring when taking an axillary temperature. When taking an axillary the temperature keep the arm close to the body so the probe tip is well covered by the arm. Axillary temperatures can vary from 35,2 to 36,8°C. From the medical point of view taking an axillary temperature is not a high-accuracy way of measuring. For receiving exact figures it is necessary to have a solid contact of tip with the body. When the tip is not in a close contact with the skin the unit beeps earlier and the result of measuring is lower.

• Taking an oral temperature.

Place the probe tip under the tongue as near as possible to the heat pocket. Keep the mouth closed when measuring the temperature. Oral temperatures can vary from 35,7 to 37,3°C.

⚠ **Attention!** Such activities as having hot or cold drinks, smoking, exercising may increase or decrease your temperature. Therefore it's important to relax with your mouth closed for about 5 minutes before taking an oral temperature.

• Taking a rectal temperature.

Apply a water-soluble lubricant to the probe cover. Gently insert the probe (around 1 cm) into the rectum. Rectal temperatures can vary from 36,2 to 37,7°C.

Average measuring time for different measuring methods

	WT-03	WT-04	WT-06
Axillary	2-3 min	50-60 sec	20 sec
Oral	1 min	40-50 sec	10-15 sec
Rectal	1 min	30-40 sec	6-10 sec

Taking an axillary temperature may take longer time because of the features described above. When the

temperature rises slightly the unit beeps.

⚠ **Attention!** When a unit beeps it doesn't mean that measuring is completed. It only means that your temperature is rising slightly. Keeping the unit for 30-60 sec after the beep is recommended.

Beeps definition:

	WT-03	WT-04	WT-06
Short	t ≥ 37,5°C	t ≥ 37,8°C	t ≥ 37,8°C
Long	t < 37,5°C	t < 37,8°C	t < 37,8°C

In order to prolong the battery life, turn OFF the unit after use pressing ON/OFF button. The unit will automatically turned OFF approximately after 10 minutes.

ⓘ **Note:** According to hygiene standards do not use for rectal measuring the same unit as used for axillary or oral measuring. Avoid having hot or cold drinks, exercising, smoking or taking a shower directly before measuring.

GENERAL CARE

- Children must use the thermometer under the supervision of adults.
- Do not walk, run or talk during temperature taking.
- Wipe the unit before and after each measurement.
- Store the unit in the protective case when it is not used.
- Do not gnaw the probe tip.
- Do not store the unit under direct sunlight, in dust or humid place. Avoid extreme temperatures.
- Dropping or subjecting your thermometer to strong shocks should be avoided.
- Do not attempt to disassemble the unit, except for battery replacing.
- If the display shows "ERROR" and the measurement doesn't start please contact a service center.

ⓘ **Note: The unit may be damaged if**

- it is used and/or stored out of determined temperature and humidity conditions;
- it was dropped or had a mechanical injury.

⚠ **Contraindications:** are not revealed.

REPLACING THE BATTERY

The battery in your digital thermometer needs to be replaced when " " or " " appears in the lower right corner of the display.

To replace the battery:

1. Turn off the screw of the battery cover (WT-06 only).
2. Take the battery cover away from the unit.
3. Replace the battery.
4. Shut the battery cover back into place.

ⓘ **Note:** Please keep the battery away from children. Do not store the battery at a high temperature.

SPECIFICATIONS

	WT-03	WT-04	WT-06
Measurement range	32,0 °C - 42,9 °C	32,0 °C - 42,9 °C	32,0 °C - 42,9 °C
Display resolution	0,1 °C	0,1 °C	0,1 °C
Accuracy	Accuracy +0,1°C for measurement range 35,5°C - 42,0°C Accuracy +0,2°C for measurement range lower than 35,5°C and higher than 42,0°C		
Auto power OFF	10 min ± 2 min	8 min 40 sec ± 2 min	10 min ± 2 min
Battery	LR41 / SR41	LR41 / SR41	LR41 / SR41
Life of battery	Min 1 year	Min 1 year	Min 1 year
Unit size	128,1×18,5×10,5 mm	133×21×11 mm	148×22×20 mm
Weight (with battery, without case)	- 10 g	- 10 g	- 14 g
Operating condition	Temperature: 10°C - 40°C; relative humidity: 30% - 85%		
Storage condition	Temperature: -10°C - 60°C; relative humidity: 25% - 90%		

THERMOMETER SET

1. Thermometer
2. Battery
3. Instructions and warranty card
4. Protective case

UTILIZATION

The unit must be utilized in accordance with current standards separately from domestic wastes. For utilizing it is necessary to contact special organizations licensed to make utilization.

MEETING STANDARDS

This unit meets EC requirements for medical equipment 93/42/EEC. Declaration of conformity.

MAINTENANCE

It is recommended the performance should be checked every 2 years or after it has been dropped. Please contact the service center.

ⓘ **Note:** You may check the unit accuracy yourself. Pour a vessel (300-500 ml) with warm water from 36 to 41°C. Put two or more thermometers including one you want to check and a reference one.

Probe tips must be at the same level and 2-3 cm in water.

After 1-2 min compare the figures. Their difference must be not more or less than 0,1°C.

WARRANTY

Warranty period is 2 years from the date of purchase. This warranty doesn't cover any damages caused by improper using, and also battery, protective cover and packaging. When a manufacturing defect is revealed during the warranty period a faulty unit would be repaired or, if repairing is impossible, replaced with another one.

The date of manufacture is indicated on the device in the serial number SN: WYYXXXXX. The first and second digits (WW) are the week of production, the third and fourth (YY) are the year of production.

ⓘ **The manufacturer may change units partially or completely if necessary, without prior notice.**

SYMBOL INFORMATION

<p>IP22 FOLLOW INSTRUCTIONS BEFORE USE HOUSING INGRESS PROTECTION RATE-IP22 (Protected from the penetration of solid bodies with dimension greater than 12.5 mm. Protected from the penetration of vertically falling water drops)</p>	<p>CE CE MARK</p>	<p>DISPOSAL FOR SEPARATE COLLECTION TYPE BF EQUIPMENT ARTICLE NUMBER LOT NUMBER OPERATING CONDITION TEMPERATURE STORAGE CONDITION TEMPERATURE</p>
---	------------------------------	---

Last revision 2021-W08

B.Well Swiss AG
Bahnhofstr. 24, 9443 Widnau, Switzerland
www.bwell-swiss.ch