EN Instructions for Use 4

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



6-4

1. INTRODUCTION

Thank you for purchasing the B.Well upper arm blood pressure monitor PRO-25.

PRODUCT DESCRIPTION, INTENDED USE AND APPLICATION AREA

The arm blood pressure monitor is intended to measure the systolic pressure and diastolic pressure, as well as the pulse rate of adult person via non-invasive oscillometric technique at medical facilities or at home. The device is not intended for use on infants and children.

The device is designed for home or clinical use.

Measurement position is on upper arm only. Intended patient population: >12 years of age, i. e. adolescents age 12 through 18 years of age, adults age >18 years of age.

- ⚠ In case of any serious incident with the device, immediately contact the manufacturer and national competent authorities.
- ⚠ Lay person or clinical professionals. Can read and understand the user manual.
- ⚠ For adult only except for whom with contraindications.
- ① 1. The patient is an intended operator, and all the functions can be safely used.
- 2. Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.

IMPORTANT ADVANTAGES OF PRO-25

- Accuracy control indicators: cuff loose, arm shake.
- Universal cuff size 22–42 cm fits almost every arm.
- Correct monitoring as recommended by the doctor – reliable blood pressure level through averaging of 3 readings.
- Memory for 2 users of 99 values each with date and time.
- Blood pressure classification scale: easy interpretation.
- Pulse arrhythmia indicator.
- Ergonomic design and soft screen protection.
- USB-C cord.

△ PRE-APPLICATION INFORMATION

Please read this manual thoroughly before using the unit. Please retain this manual for future reference. For specific information about your blood pressure, please CONSULT YOUR DOCTOR. To avoid risk and damage follow all warning precautions. Operate unit only as intended. Read all instructions prior to use.

The use of third-party components (cuffs, air tubes and connectors) may negatively affect the operation of the device, lead to failure of the device algorithm and cause harm to the user.


Batteries, USB cables and adapters, if required when using this device, must comply with the requirements specified in the instructions for use.

Use of unsuitable power supply elements (batteries, USB cables) may cause electromagnetic interference and impair operation of the device.

2. CONTENTS AND DISPLAY INDICATORS

Watch screen, device and cuff pics (see pictures on pages 2 and 3).

LCD SCREEN 2-1

- ① User 1
- ② User 2
- ⋈ Inflation icon
- ⋈ Deflation icon
- ⋈ Memory icon • Memory number

-  Cuff loose indicator
-  Arm shake indicator
-  Heartbeat icon
-  Pulse arrhythmia indicator
-  Average blood pressure icon
-  Battery icon
-  ESH* Blood pressure classification indicator

28:28 Date and time

DEVICE 2-2 2-3

- 2-2-1** Display
- 2-2-2** User selection/Date/Time/Setting button
- 2-2-3** START/STOP button
- 2-2-4** Memory button
- 2-2-5** USB-C jack
- 2-3-1** Air jack

CUFF 2-4

- 2-4-1** Arm cuff
- 2-4-2** Air plug
- 2-4-3** Air tube

3. CONTRAINDICATION

- ⚠ Do not use this device if the patient's condition meets the following contraindications, to avoid inaccurate measurements or injuries.
- ⚠ The device is not suitable for use on patients with implanted, electrical devices, such as cardiac pacemakers, and defibrillators.
- ⚠ Avoid taking measurement on the arm on the side of a mastectomy or lymph node clearance.
- ⚠ The device measures blood pressure using a pressurized cuff. If the measuring limb suffers from injuries (for example open wounds) or under conditions or treatments (for example intravenous drip) making it unsuitable for surface contact or pressurization, do not use the device, to avoid worsening of the injuries or conditions.

* European Society of Hypertension.

- ⚠ Avoid taking measurements of patients with conditions, diseases, and susceptible to environment conditions that lead to uncontrollable motions (e.g. trembling or shivering) and inability to communicate clearly (for example children and unconscious patients).
- ⚠ The device uses oscillometric method to determine blood pressure. The arm being measure should have normal perfusion. The device is not intended to be used on a limb with restricted or impaired blood circulation. If you suffer with perfusion or blood disorders, consult your doctor before using the device.

4. ABNORMAL REACTION

Not detected.

5. PRECAUTIONS AND WARNINGS

5.1. Precautions before and during use

- ⚠ Before use, please read this instruction manual carefully.
- ⚠ DO NOT use this monitor in a moving vehicle such as in a car.
- ⚠ Motion, trembling, shivering may affect the measurement reading.
- ⚠ DO NOT use this monitor for infants, toddlers, children or persons who cannot express themselves.
- ⚠ DO NOT use this monitor for any purpose other than measuring blood pressure and pulse rate.
- ⚠ DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.
- ⚠ DO NOT use in a location where there is moisture or a risk of water splashing this monitor. This may damage this monitor.
- ⚠ DO NOT drop or subject this monitor to strong shocks or vibrations.
- ⚠ Do not use or store the monitor outside the manufacturer's specified conditions (extremely high or low temperatures and humidity), as this may affect the performance or cause inaccurate measurements.
- ⚠ When the performance changes (such as: inaccurate measurement or abnormal display),

please stop using it immediately and contact the sales service personnel in time.

- ⚠ Do not perform measurements more frequently than necessary. Due to the interference of blood flow, some bruising may occur.

5.2. Medical precautions

- ⚠ Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
- ⚠ DO NOT take medicine based on readings from the device. Contact your physician for specific information about your blood pressure. The patient should not self-diagnose or self-medicate per measured results. Kindly adhere to the instructions of your physician or health provider.
- ⚠ DO NOT use the device while you are on an intravenous drip or blood transfusion.
- ⚠ Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular therapy, or an arteriovenous (AV) shunt.
- ⚠ Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.
- ⚠ For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
- ⚠ Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion, diabetes, pregnancy, pre-eclampsia or renal disease.
- ⚠ Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.
- ⚠ People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.

- ⚠ The cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
- ⚠ Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders, because the cuff inflation can cause bruising.

5.3. General precautions

- ⚠ When the ambient temperature is less than 5°C, please take the device to the place where the ambient temperature is between 5°C–40°C at least 1 hour. When the ambient temperature is higher than 40°C, please take the device to the place where the ambient temperature is between 5°C–40°C at least 2 hours.
- ⚠ DO NOT use this monitor in areas containing high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, computerized tomography (CT) scanners. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.
- ⚠ Maintenance should be done by the manufacturer as suggested.

5.4. Warning / Safety notice

- ⚠ Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self treat these symptoms without consulting your physician first.
- ⚠ If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
- ⚠ Consult with your physician before using this monitor on an arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow which could result in injury.
- ⚠ Consult with your physician before using this monitor if you have had a mastectomy or lymph node clearance.
- ⚠ Product is designed for its intended use only. Do not misuse in any way.
- ⚠ Do not use the monitoring ME EQUIPMENT on the same limb simultaneously. This could temporarily cause loss of function or an inaccurate measurement.

- ⚠ Too frequent measurements can cause injury to the PATIENT due to blood flow interference.
- ⚠ DO NOT attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
- ⚠ The cuff should not be applied over a wound as this can cause further injury.
- ⚠ Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
- ⚠ A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.
- ⚠ Please check whether the operation of the Arm blood pressure monitor leads to prolonged impairment of patient's blood circulation by observing the limb concerned.
- ⚠ Use only the approved arm cuff for this unit. Use of other arm cuffs may cause incorrect measurement results.
- ⚠ The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity range.
- ⚠ Do not use the device near strong electrical or electromagnetic fields generated by cell phones or other devices, they may cause incorrect readings and interference or become interference source to the device. Do not use the device during patient transport outside healthcare facility for interference source existing as well.
- ⚠ No maintenance or servicing when using.
- ⚠ Do not mix new and old batteries simultaneously.
- ⚠ Do not mix battery types. Do not insert the batteries with their polarities incorrectly aligned. Long-life alkaline batteries are recommended.
- ⚠ Only use a recommended AC adapter double-insulated complying with en 60601-1 and en 60601-1-2. An unauthorized adapter may cause fire and electric shock. It is recommended to use B-Well adapter: AD-75.
- ⚠ Please place the blood pressure monitor out of the reach of children.
- ⚠ Please use component (e.g. cuff) provided by manufacturer. Otherwise, the measurement accuracy will be affected.
- ⚠ No modification of this equipment is allowed.

- ⚠ To avoid strangulation, please keep the air tube and USB-C charging cable away from the infants, toddlers and children.
- ⚠ Do not leave the small parts where children can reach them. Children may swallow them. If a child accidentally swallows them, battery cover, please contact a doctor immediately.
- ⚠ The cuff complies with the requirements of ISO 10993-5, ISO 10993-10, ISO 10993-23. But few sensitive people may have allergies.
- ⚠ DO NOT use this monitor on an injured arm or an arm under medical treatment.

6. SETUP AND OPERATING PROCEDURES

6.1. Battery loading and replacement

After unpacking the device install the batteries. The battery compartment is located on the underside of the case.

Battery installation:

1. Slide battery cover off as indicated by arrow on the cover.
2. Install 4 new AAA alkaline batteries according to polarity.
3. Close battery cover.

6-1

- ⚠ Do not leave discharged batteries inside the device. Also remove the batteries if the device will not be used for a month or more to avoid possible damage to the device in case of battery leakage.
- ⚠ Avoid contact of electrolyte with eyes. If electrolyte gets into eyes, rinse immediately with plenty of clean water and consult a doctor.
- ⚠ If the product won't be not used for a long time (over 3 months), please take out the batteries. Dispose batteries properly; observe local laws and regulations.
- ⓘ The adapter access function is suitable for temporary use when you don't have suitable AAA batteries around. The adapter should comply with the requirement of IEC 60601-1 standard and the specifications must meet the requirements:
INPUT: AC 100~240 V, 50/60 Hz,
OUTPUT: DC 5 V, 1.0 A.
Other AC adapter may vary in output voltage and

polarities and may represent a risk on your life and damaging the device.

Low battery indicator

If the batteries are low, the corresponding symbol  will appear on the screen. The appliance will then not switch on. Replace all batteries with new ones.

6.2. Date and time setting

In standby mode, press the  button for about 3 seconds to enter the date setting, and "year" will flash. Press the  button to adjust to the desired year, and then press the  button to confirm the selection. When the year is set, it will automatically enter the month setting. At this time, the month icon will flash. You can switch to the desired value by pressing the  button, and then press the  button to confirm the selection. Follow the same step to set date, hour and minute.

6-2

6.3. Connecting components

Connecting the AC adapter / cable

The device has a USB-C connector for connecting the AC power adapter / cord, which is located at the right side panel of the device.

1. Plug the USB-C adapter / cord into a 100–240 V, 50–60 Hz AC mains supply.
2. Connect the plug of the AC adapter / cord to the appropriate connector.

The battery power of the unit is automatically switched off.

- ⓘ If the AC adapter or USB-C cord is not included, it can be purchased separately. Use only the AC adapter / cord recommended by B.Well. If any other adapter is used, the service centre will not be liable for warranty service.

- ⚠ Do not plug or unplug the power cord into the electrical outlet with wet hands.
- ⚠ Do not overload power outlets. Plug the device into the appropriate voltage outlet.
- ⚠ If the AC adapter is abnormal, please change the adapter.
- ⚠ Do not use any other type of AC adapter as it may harm the monitor.
- ⚠ Please use the adapter with characteristics:
INPUT: 100–240 V–50/60 Hz, 0.2 A;
OUTPUT: 5.0 V, 1.0 A, 5.0 W.

Connecting the arm cuff to the monitor

Firmly insert air tube plug into opening located on the left side of monitor unit. Make certain that the plug is completely inserted in order to prevent air leakage during use. With sticky nylon section facing outward, insert end of cuff underneath metal ring of cuff.

- ⚠ Avoid compression or restriction of the connection tubing during measurement, which may cause inflation error, or harmful injury due to continuous cuff pressure.

6.4. Applying the cuff

- a) Thread the end of the cuff through the metal ring so that the entry point of the air tube into the cuff is on the outside.
- b) Slide your arm through the cuff so that the air tube is positioned along your arm towards your palm.
- c) Fasten cuff about 1–2 cm (0.4–0.8") above the elbow joint.
- d) Tighten the cuff firmly but not too tightly, taking into account the taper of the arm, by pulling the free end.
- e) The cuff must fit tightly around the arm, otherwise the measurement result will be incorrect. The cuff must be at the level of the heart. It is not recommended to wear the cuff over clothing.
- f) Place your hand with the palm of your hand facing upwards. Make sure the pump tube is not twisted, compressed or restricted.

6-3

6-4

- ⚠ The Artery mark on the cuff indicates the location of the artery on the left arm, as this is the most common arm for measuring blood pressure.

6.5. Carrying out a measurement

Before measuring

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position while keeping your arm level with your heart. Do not cross your legs.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.

6. Wait 5 minutes or longer before re-testing.
 7. Try to measure your blood pressure at the same time each day for consistency.
 8. Test comparisons should only be made when monitor is used on the same arm, in the same position, and at the same time of day.
 9. This blood pressure monitor is not recommended for people with severe arrhythmia.
 10. Do not use this blood pressure monitor if the device is damaged.
- ⚠ The original cuff should be used. Any replacement of the original component not provided by the manufacturer may cause an incorrect measurement result.

Measurement in sitting position

1. Sit down and place your feet flat on the floor, don't cross your legs or bend them under you. Lean your back straight against the back of the chair.
 2. Place palm upside in front of you on a flat surface such as table.
 3. The middle of the cuff should be at the level of the heart.
- ⚠ To avoid errors in the measurement, it is important to remain immobile during the measurement and in silence.

Measuring in the prone position

1. Lie on your back.
2. Place your left arm straight along your side with your palm upside.
3. The cuff should be placed at the same level as your heart.

Common sources of errors or incorrect pressure values

1. Movement during measurement.
2. Incorrect positioning of the cuff in relation to the heart.
3. The cuff does not fit you in size.
4. Loosely tightened cuff.

Power on

Press the  button to turn the unit on. All the existing symbols will appear on LCD screen for one second as unit performs a quick diagnosis.

Measuring

Press the  button. After all icons are turned on, the monitor will start inflating for measurement and display . Check the measured values after the measurement finished.

- ① If you feel uncomfortable during the measurement press the  button immediately to stop measure.
- ① When the air pressure is filled to a certain value, the value on the display screen will slowly drop at a certain speed and  symbol will flash. After the measurement is completed, the blood pressure values and pulse measurement will be displayed on the screen. An indicator representing the current measurement will appear next to the corresponding ESH Blood pressure classification scale.

6-4

Memory function

Each measured value is stored automatically under the appropriate user group. This device can store up to 99 sets of measurements for each of 2 users. Once the memory log is full, old values will be refreshed with new ones.

Average value of last 3 results setting

In the power-off mode, press the  button once and the device will display the average value of the blood pressure measurements of the last 2 or 3 times, and the **AVG** indicator will appear. Press the  button again, and the latest measured value will be displayed. Press the  button again and the rest measured values will be displayed one by one.

Delete memory

Press the  button to power off the device. Press the  button to select the user group whose measured values need to be deleted. Press the  button again to power off the device. Hold down the  button for about 5 seconds. When display flashes the  icon, press the  again to clear the memory and turn off the device

Memory group settings

Press the  button to power on the device and press button  to select the memory group and switch between user 1 or user 2.

Pulse arrhythmia indicator

When the monitor detects an irregular rhythm during the measurements, the Pulse arrhythmia indicator  will appear on the display with the measurement values. Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure.

- ① Consult your physician if the Pulse arrhythmia indicator frequently appears with your test results.

TAKING CARE OF MEASUREMENT

ACCURACY:

Arm shake Indicator

The  icon flashes when you move your body or shake your arm during the measurements, which may cause incorrect measurement results. Please, adjust your posture and measure again.

Cuff loose Indicator

The  icon is always displayed on the screen when the cuff is wrapped correctly. When the cuff is too loose, the  icon will always flash to remind you. If the  icon is flashing, please, press the  button to stop the measurement.

Power off

Press the  button to turn off the arm blood pressure monitor. The monitor automatically turns off after 1 minute of inactivity.

6.6. Blood Pressure Classification scale: easy interpretation of measurement results

The following guideline for assessing high blood pressure (without regard to age or gender) has been established by the European Society of Hypertension (ESH). Please note that other factors (e.g. diabetes, obesity, smoking, etc.) need to be taken into consideration.

Range	SYS mmHg	DIA mmHg	Color indicator	Measures
Optimal	< 120	< 80	Green	Self-check
Normal	120–129	80–84	Green	Self-check

High normal	130–139	85–89	Green	Consult your doctor
Grade 1: mild hypertension	140–159	90–99	Yellow	Consult your doctor
Grade 2: moderate hypertension	160–179	100–109	Orange	Consult your doctor immediately
Grade 3: severe hypertension	≥ 180	≥ 110	Red	Urgently seek medical advice
ESH Definition and Classification of Blood Pressure Levels				
SYS – systolic blood pressure, DIA – diastolic blood pressure				

- ⚠ The blood pressure level indication is a reference and does not imply independent emergency measures.
- ⚠ Consult with your physician for accurate assessment, and never change your treatment by yourself.
- The colored scale shows the range within the measured blood pressure value is located.
- Depending on the height of the indicator, the readout value is either within the normal (green), borderline (yellow and orange) or hazardous (red) range.
- If the systole and diastole values fall into two different categories (e.g., systole in the High normal category and diastole in the Normal category), the display always shows the higher category.

7. TROUBLESHOOTING

- If an error occurs during measurement, the measurement process will be interrupted, the instrument will sound an error message and prompt you to proceed.
An error message (e.g. Er 1) appears on the display.

Troubleshooting (1)

LCD SYMBOLS	POSSIBLE CAUSE	SOLUTION
Er U	The inflation can not reach 30 mmHg in 12 seconds	Please fasten the cuff correctly and measure again
Er H	The inflation reaches 295 mmHg	
Er 1	The pulse rate is not detected correctly	

Er 2	Too much disturbance (move, talk, or magnetic disturbance during a measurement)	Please measure again. Avoid speaking or moving body parts while measuring
Er 3	The measurement result is abnormal	Please fasten the cuff correctly and measure again. If you cannot solve the problem, please contact the manufacturer
Er 23	SYS value is lower than 57 mmHg	
Er 24	SYS value is higher than 255 mmHg	
Er 25	DIA value is lower than 25 mmHg	
Er 26	DIA value is higher than 195 mmHg	

- If errors occur during measurement, check the following points and take appropriate action:

Troubleshooting (2)

PROBLEM	POSSIBLE FAULTY	SOLUTION
Failure to power on	Whether the power is insufficient	Replace the batteries or insert the USB-C charging cable for power supply
	Whether the positive and negative poles of the battery are installed reversely	Install the batteries correctly
No pressurizing	Whether the air tube plug is inserted tightly	Insert the air tube plug firmly into the jack
	Whether the air tube is broken or leaked	Please contact the dealer to replace with a new cuff
Unable to measure due to the display error	Whether the arm is moved when pressurization	Keep your arm and body still
	Whether you talk during the measurement	Keep quiet while measuring the blood pressure
Air leakage of the cuff	Whether the cuff is twined too loose	Please tighten the cuff
	The airbag of the cuff is ripped	Please contact the dealer to replace with a new cuff

- △ If the blood pressure still cannot be measured after trying the stated solutions, please contact the dealer. Do NOT attempt to disassemble the device by yourself.

8. MAINTENANCE

8.1. Maintenance (general)

- △ Avoid dropping, slamming, or throwing the unit.
- △ Avoid high temperature and direct sunlight. Do not immerse the monitor in water as this will result in damage to the monitor.
- △ Do not attempt to disassemble this monitor.
- △ Do not allow water or other liquids to flow into the device. The arm pressure monitor should not no longer be reused when liquid enter and damage the device and cuff.
- ① Do not bend or crease the air tube excessively.
- ① Do not store the monitor or its components:
 - if the monitor or its parts is wet;
 - in locations with extreme temperatures, humidity, direct sunlight, dust, or corrosive gases;
 - in areas with a high risk of vibrations or shocks.

8.2. Maintenance (operation)

- Remove batteries when not in operation for an extended period of time.
- It is recommended the performance should be checked every 2 years of regular use or after dropping/damaging. Please contact an authorized service center.
- △ Do not kink the connecting tubing, as the resulting continuous cuff pressure can cause interference with blood flow and harmful injury to the patient.
- ① The bottom edge of the arm cuff should be 0.79–1.18 inches (2 to 3 cm) above the inside elbow. The air tube is on the inside of your arm and aligned with your middle finger.
- ① Repeated measurement will result in blood congestion in the arm, which will affect the measurement result.
- ① Consult your doctor if unexpected readings are obtained.

8.3. Maintenance (cleaning)

- △ When cleaning the unit, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent. Do not use solvents.
- △ **Cuff cleaning:** it is recommended that the cuff be disinfected before use for different users. Wipe the inside of the cuff (skin contact side) with a soft cloth wrung out after soaking in ethyl

alcohol (75–90%) or 3% hydrogen peroxide solution, then air dry the cuff.

- ⚠ Do not disinfect through methods like high-temperature steam or ultraviolet radiation. These might damage the device and reduce its service life.
- ⚠ The cuff cover can be washed at 30°C. Remove the bladder beforehand through the special technological hole in the cuff cover.
- ⚠ The bladder must not be washed! The cuff cover must not be ironed or dried in the tumble dryer! The technological hole must not be sewn up!
- ⚠ Do not use any corrosive cleaning agent. When cleaning, be careful not to immerse any part of the monitor to avoid liquid flow into the instrument.

The device, batteries and cuffs must be disposed of in accordance with local regulations at the end of their useful life.

9. SPECIFICATIONS

Product description	Digital upper arm blood pressure monitor
Model	PRO-25
Display	LCD size: 44.5 mm × 58.5 mm (1.75" × 2.30")
Measuring method	Oscillometric method, automatic air inflation and measurement
Classification	Internal powered equipment, Type BF cuff is the applied part 
Unit dimensions	Approx. 118 mm × 98 mm × 62.5 mm (4.65" × 3.86" × 2.46")
Cuff circumference	22–42 cm (8.6"–16.5")
Unit weight	Approx. 220 g (excluding batteries)
Memory	2 × 99 memories with date and time
Power source	Batteries: 4 × 1.5 V  AAA DC: 5.0 V  1.0 A Medical AC adapter: INPUT: 100–240V-50/60 Hz, 0.2 A; OUTPUT: 5.0 V, 1.0 A, 5.0 W
Measurement range	Cuff pressure: 0–295 mmHg SYS: 57–255 mmHg DIA: 25–195 mmHg Pulse rate: 40–199 beats/minute
Accuracy	Pressure: ±3 mmHg Pulse rate: ±5%

Operating conditions	Temperature: 5°C–40°C (41°F–104°F) Humidity: 15%–90% RH (non-condensing) Atmospheric pressure: 70KPa–106KPa
Transportation & storage conditions	Temperature: –20°C–55°C (–4°F–131°F) Humidity: 10%–93%RH Atmospheric pressure: 70KPa–106KPa
Battery life	Approx. 2 months at 3 tests per day

① These specifications are subject to change without notice.

10. PACKAGE CONTENTS

- Upper arm blood pressure monitor.
- Universal M–L cuff (upper arm circumference 22–42 cm), article C-ML-2242-ST.
- AAA batteries – 4 pcs.
- Mains USB adapter, article AD-75 (depending on the set).
- USB-C power supply cord (depending on the set).
- Storage bag.
- Gift box.
- Instruction manual.

11. SYMBOL INFORMATION

	FOLLOW INSTRUCTIONS BEFORE USE
	WARNING
	INFORMATION
	MANUFACTURER'S NAME
	ARTICLE NUMBER
	SERIAL NUMBER
	CE MARK WITH NUMBER OF NOTIFIED BODY
	BF TYPE DEVICE
	DISPOSAL FOR SEPARATE COLLECTION
	RECYCLABLE

IP21	HOUSING INGRESS PROTECTION RATE: IP21
	KEEP DRY
	KEEP AWAY FROM SUNLIGHT
5 °C  40 °C OPERATION CONDITION	OPERATING CONDITION, TEMPERATURE
-20 °C  55 °C STORAGE CONDITION	STORAGE CONDITION, TEMPERATURE
 90 % 15 %	OPERATING CONDITION, HUMIDITY
 93 % 10 %	STORAGE CONDITION, HUMIDITY
 106 kPa 70 kPa	OPERATING CONDITION, ATMOSPHERIC PRESSURE
MD	MEDICAL DEVICE
UDI	UNIQUE DEVICE IDENTIFICATION
EC REP	EC REPRESENTATIVE

12. WARRANTY INFORMATION

Warranty period is 3 years from the date of purchase for monitor and 1 year for cuff and adapter.

This warranty doesn't cover any damages caused by improper using, and also battery, and packaging. When a manufacturing defect is revealed during the warranty period a faulty unit would be repaired or, if repairing is impossible, replaced with another one. The warranty does not cover components and consumables subject to wear and batteries, bags, and package of the item.

Manufacturing date is in a serial number: WWWWXXXXX. WW – week, YY – year of production.

XXXXX is the consecutive number in the batch.

The manufacturer may change units partially or completely if necessary, without prior notice.

13. ELECTROMAGNETIC COMPATIBILITY INFORMATION

⚠ Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility.

■ The essential performance:

1. Limits of the error of the cuff pressure indication.
2. Reproducibility of the blood pressure DETERMINATION.
3. Alarm.

When electromagnetic interference affects the above performance, please stop using the device.

■ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

■ Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

■ It is suggested that the blood pressure monitor be kept at least 30 cm away from other wireless devices, such as WLAN unit, microwave oven, etc. It can't be used near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.

14. ADDITIONAL NOTES

- ① When measuring with the cuff, if the tester feels seriously uncomfortable, press the button of the blood pressure monitor to deflate the cuff, or remove the cuff directly from the arm.
- ① Do not stand while testing. Sit in a relaxed position while keeping your arm level with your heart.
- ① Avoid speaking or moving body parts while testing.
- ① When cleaning the unit, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent.

- ① Cuff cleaning: Do not soak cuff in water! Apply a small amount of rubbing alcohol to a soft cloth to clean cuff's surface. Use a damp cloth (water-based) to wipe clean. Allow cuff to dry naturally at room temperature. The cuff must be cleaned and disinfected before use between different users.
- ① It is recommended the performance should be checked every 2 years.
- ① The equipment can be cleaned by lay operator according to the cleaning procedures in the instructions.
- ① The monitor, the batteries and the cuff, must be disposed of according to local regulations at the end of their usage.

Last revision 2025-W14

B. Well Swiss AG,

Bahnhofstrasse 24, 9443 Widnau, Switzerland

www.bwell-swiss.ch

B.Well 
caring for everyone

B.Well Swiss AG,

Bahnhofstrasse 24, 9443 Widnau, Switzerland

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B. Well Swiss AG
Bahnhofstrasse 24, 9443 Widnau,
Switzerland
www.bwell-swiss.ch

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